



Empore™

Extraction Disk Plate Accessories

PRODUCT _____
Information

Tape Pad

General Information

3M™ Empore™ Tape Pads are intended for sealing unused portions of a 96-well solid phase extraction plate during vacuum processing to enable vacuum to be uniformly applied to each well being used. Use a clean sheet each time for maximum sealing ability. **Note:** these tape sheets are not intended for sealing collection plates prior to sample injection because the adhesive may clog the injection needle.

Product Information

The Empore tape pad (Catalog # 660) is a pre-cut tape pad of high performance tape sheets. Tape sheet construction is a polypropylene film coated with a synthetic rubber adhesive. Each sheet is coupled with a liner tab for easy handling and dispensing.

Product Characteristics

Sheets/Pad:


25 each

Pads/Pack:

10 each

Size:

5.5" x 3.25"



Application Technique

Sheets of Empore tape pad should be removed one at a time by grasping the end with the liner and pulling the sheet free from the pad. Place the far end of the sheet in the desired location and allow the sheet to adhere from that end back to the tabbed end. Remove the tab and press the whole sheet flat. Application of the tape sheets is easiest at room temperature. At colder temperatures, approaching 32°F (0°C), the adhesive becomes more firm and may not adhere as well.

Warranty Information: All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M WARRANTS ONLY THAT ITS PRODUCTS WILL MEET 3M SPECIFICATIONS AT THE TIME OF SHIPMENT. THE FOREGOING WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND FREEDOM FROM NON-INFRINGEMENT. **Empore™ Sample Preparation Products are intended for solid phase extraction during scientific research only. These products are not intended or warranted for use in medical devices or in assessment and treatment of clinical patients.**

Limitation of Remedies: If Products are proven not to meet 3M specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: THE REMEDIES PROVIDED HEREIN ARE EXCLUSIVE REMEDIES AGAINST 3M FOR ANY ALLEGED OR ACTUAL NONCONFORMANCE TO SPECIFICATIONS OR DEFECT OR OTHER FAILURE IN PRODUCTS. UNDER NO CIRCUMSTANCES IS 3M LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING LOST PROFITS) IN ANY WAY RELATED TO THE PRODUCT OR SUPPLY OF PRODUCT UNDER ANY THEORY OF LAW INCLUDING, BUT NOT LIMITED TO CONTRACT, NEGLIGENCE AND STRICT LIABILITY.

3M Bioanalytical Technologies Project

3M Center, Building 270-02A-08
St. Paul, MN 55144-1000
1-888-509-5330
Fax: 1-651-736-4882
Website: www.3M.com/empore