

3M™ Cavilon™ NSBF*



No sting barrier film

Product description

- ▶ A sting-free, alcohol-free liquid barrier film that dries quickly to form a breathable, transparent coating on the skin. Designed to protect intact or damaged skin from urine, feces, other body fluids, tape trauma, and friction.

Features and benefits

- ▶ **Alcohol free:**
will not sting even when applied to damaged or denuded skin.
- ▶ **Provides superior protection from adhesive tape trauma:** it reduces skin tears and complications related to tape use.
- ▶ **Reduces time needed for wound care:**
under normal use, reapply once every 72 hours.



Indications

Protecting intact or damaged skin from body fluids, adhesives and friction:

Suggested applications include:

- ▶ Anal wounds (tumor removal).
- ▶ Perineum wounds (hernia).
- ▶ Vulvar wounds (episiotomy).
- ▶ Stomy wounds (urethrostomy, trans-cutaneous gastric tubing).
- ▶ Chops wounds
- ▶ Mud fever (in horses).
- ▶ Lowering of auditory duct.
- ▶ Dressing for very small pets.

* No Sting Barrier Film.



Laboratoires 3M Santé
Département Vétérinaire
Boulevard de l'Oise
95029 Cergy Pontoise Cedex
France
Tél. : 33.(0)1.30.31.82.82
Fax : 33.(0)1.30.31.83.83
www.3MSante.fr

3M



Wound management

Wound management

3M™ Clear Custom Collar



Product description

- ▶ Clear collar, with padded neck edge, that can be custom fit to each individual animal's neck, thanks to the 3M Scotchmate™ Hook and Loop Fastener system.

Features and benefits

- ▶ **Clear:** greater visibility for the animal means more comfort and less disorientation. The owner is satisfied.
- ▶ **Customized fit:** custom fit to each individual animal, no need to fasten to the animal's collar.
- ▶ **Easy and quick application:** collar may be removed and re-applied easily, even by the owner.
- ▶ **Comfortable:** thanks to the padded neck edge, no risk of wound or loss of hair.
- ▶ **Light, flexible and resistant.**

3M™ Vetbond™

Tissue adhesive



Product description

- ▶ N-butyl cyanoacrylate adhesive polymerizes in seconds after contact with tissue and body fluids.

Features and benefits

- ▶ **Quick wound repair:** 3M Vetbond tissue adhesive may be used in a variety of vet applications to bind wound edges together, hence reducing the need for classical sutures and analgesia in minor wound repair.
- ▶ **It reduces tissue trauma** and increases the wound appearance through decreasing scar marks.
- ▶ **As healing occurs, 3M Vetbond is sloughed off naturally.**



Indications

- ▶ Suture of minor wound (less than 3 cm long). If the wound is longer, support with staples.
- ▶ Suture of small surgical incisions (pins removal in cats) or as complement of subcutaneous continuous sutures (cats ovariectomy, minor wounds under bandaging tapes during articular surgeries).
- ▶ Corneal ulceration smaller than half the stroma width.
- ▶ Oral surgery (teeth removal).
- ▶ Suture of minor wound in very small pets (rodents, exotics...).

3M™ Precise™

Disposable skin stapler



Product description

- ▶ 3M Precise range of disposable skin staplers and staple removers provides the highest degrees of performance, versatility and reliability.
- ▶ 3M staplers are convenient, easy to use and disposable.

2 different staplers are available:

- ▶ **3M Vista skin stapler:** unique angled stapler head allows a clear view of the wound area from your normal vantage point to ensure precise staple implant. This stapler dispenses box-style staples.
- ▶ **3M DS skin stapler:** ergonomically and functionally designed. Versatile (may be used with or without handle). It dispenses arcuate, suture-style staples.

Features and benefits

- ▶ **Safety:** the rigid design holds wound edges securely by compressing the tissue between the legs of the staple. The stainless steel construction reduces inflammatory reactions and infection risks.
- ▶ **Short trigger reduces hand fatigue and allows quick and easy use.**
- ▶ **No need for analgesia:** Both styles of staplers offer speed of closure, elimination of cross-hatching and ease of removal with less wound and patient trauma.

Indications

- ▶ Skin closure (every simple, recent, clean wound, and after skin preparation).

Contraindications

- ▶ Wound with torn edges.
- ▶ Possible functional discomfort locations between the fingers or support areas.
- ▶ Bites or every deep wound (need for surgery).

3M™ Medipore™

Soft Cloth surgical tape



Product description

- ▶ 3M Medipore soft cloth tape offers excellent adhesion yet is gentle to the skin and conformable. It is suitable with a variety of applications, mainly on areas where movement is expected and on sensitive skins.

2 versions:

- ▶ **Pre-cut every 5 cm, no paper liner:**
easy to use, easy-tear perforated rolls. No need for scissors. Easy to handle, no curling or self adhesion. No waste.
- ▶ **With paper liner, non pre-cut.**



Features and benefits

- ▶ **Quality cloth backing with cross and diagonal stretch,** it conforms to body contours and to movements.
- ▶ **Superior adhesivity.**
- ▶ **Hypoallergenic:** gentle to the skin to reduce risk of skin irritation during removal. Convenient for fragile skin pets such as rodents and exotics.
- ▶ **Resistant to licking.**

Indications

- ▶ Anchoring post-surgical dressings.
- ▶ Large size dressings.
- ▶ Securing dressings on areas where swelling or movement is expected, or on fragile skins.
- ▶ To secure or cover devices, eg IV catheter sites (peripheral and central lines).



3M™ Medipore™ +Pad

Soft cloth adhesive wound dressing



Product description

- ▶ An all-in-one sterile dressing that features a breathable soft cloth tape backing with a unique non-adherent pad. It is suitable for post-surgical dressings and draining wounds.

Features and benefits

- ▶ **Highly conformable dressing** conforms to body contours, and flexes with movement or swelling. Dressing is porous and breathable to reduce risk of skin maceration.
- ▶ “Round” corners improve **adhesion** along time.
- ▶ **Hypoallergenic adhesive.**
- ▶ **Absorbent pad** protects wound, absorbs draining, and promotes wound healing, with less pain and trauma during removal.
- ▶ **Full range of sizes** is available to fit every application.

Indications

- ▶ Light to moderately draining wounds.
- ▶ Cat ovariectomy.
- ▶ Dog gelding.

3M™ Durapore™

Silk like Surgical tape



Product description

- ▶ A "silk like" tape with a hypoallergenic adhesive that offers high strength and adhesion for general purpose.

Features and benefits

- ▶ **High mechanical resistance backing.**
- ▶ **Hypoallergenic adhesive:** no residue when removed, yet strong adhesion.
- ▶ Convenient **bi-directional tear** (no need for scissors).
- ▶ Backing is **porous and breathable** for better exchange of moisture vapour and oxygen, and reduced risk of skin maceration.

Indications

- ▶ Securing variety of dressings.
- ▶ Adhesive and non-elastic support.
- ▶ Securing medical devices (peripheral catheters, drains, tubing...).
- ▶ Secural of small pets during surgeries



3M™ Micropore™

Microporous surgical tape



Product description

- ▶ An economical, general-purpose paper surgical tape that is gentle to the skin yet adheres well. It is today the benchmark microporous surgical tape for wound management.

Features and benefits

- ▶ **Adheres well:** holds well on damp skin.
- ▶ **Painless removal:** no adhesive residue on the skin
- ▶ **Porous and breathable :** for better exchange of moisture vapour and oxygen, and reduced risk of skin maceration.
- ▶ **Hypoallergenic :** gentle to the most sensitive skins.
- ▶ **Radioluscent.**



Indications

- ▶ Securing dressings:
 - when repeated taping is needed.
 - taping fragile, at-risk skin.
- ▶ Catheter secural, removed when the animal wakes up after surgery.



3M™ Microfoam™

Elastic foam surgical tape



Product description

- ▶ Elastic foam tape designed to stretch and conform while providing secure adhesion. Excellent for securing dressings and compression applications in challenging areas.

Features and benefits

- ▶ **Elastic and compressive:** secures compression dressings where movement is likely, conforms to body contours, helps to control œdema.
- ▶ **Foam backing:** allows some exchange of moisture vapour and oxygen, and protects from fluids.
- ▶ **Hypoallergenic:** gentle to the skin.

Indications

- ▶ Compression dressings (orthopedic, digestive surgical dressings).
- ▶ Post-surgical bandaging (eventration, mastectomy).
- ▶ Protection of bony prominences (under casts).
- ▶ Catheter secural.
- ▶ Articular bandaging (sprains or ligament surgery).



3M™ Transpore™

Plastic surgical tape



Product description

- ▶ A transparent, perforated plastic tape, that easily conforms to body contours.

Features and benefits

- ▶ **Microperforated:** easy straight, bi-directional tear, without scissors. Porous, it allows the skin to breathe.
- ▶ **Hypoallergenic:** good adhesive to skin yet painless when removed.
- ▶ **Transparent:** the application site remains visible, and medical devices may be easily monitored.
- ▶ **Elastic:** adheres well on sites where movement is likely.



Indications

- ▶ Securing tubing and devices (catheters, IV tubing...).
- ▶ Securing dressings.



3M™ CAT*

Cloth adhesive tape

Product description

- ▶ In addition to being natural rubber latex free and hypoallergenic, 3M Cloth Adhesive Tape has excellent adhesion for secure placement.

Product benefits

- ▶ **High adhesion:** initial, in moist and dry conditions, and along time. When removed, it leaves few residue on skin.
- ▶ **Easy to use:** easy to tear even without a pair of scissors. You can easily write on it.
- ▶ **Strong backing for durability.**

Indications

- ▶ Securing tubing, splints, dressings.
- ▶ Positioning very small pets on the table before surgery.
- ▶ Any other indications requesting high durability.

