

# Surgical Clippers and Prepping Solutions

## 3M™ Surgical Clippers by Remington and Clipper Blade Assembly

### Description

A compact and lightweight rechargeable surgical clipper with a pivoting head for removing head and body hair from patients in preparation for surgery. Uses disposable, single-use blades.

### Product Benefits

- April 1999 CDC Guidelines recommends that if hair removal is required, the use of a surgical clipper immediately before an operation is one means to help reduce the risk of surgical site infections.
- Removes all types of head and body hair, wet or dry, to allow for standardization on one product for all uses.
- Pivoting head rotates easily to allow flexibility of clipping with either a pushing or a pulling (raking) motion.
- Single-use blades reduce the risk of cross-contamination and are easy to snap on for use and pop off for disposal.
- Rechargeable clipper provides approximately 60 minutes of run time, allowing multiple patients to be clipped from a single 8-hour recharge.
- Unique design is easy to clean and disinfect.



### Suggested Applications

- Removing head and body hair from patients in preparation for procedures where hair removal is desired

	Product Number	Product Information	Items/Box	Billing Unit
<b>Packaged Sterile</b>				
	9661	Surgical Clipper	1	Case
	9662	Surgical Clipper Charger	1	Case
	9660	Surgical Clipper Single-use Blades	50	Case
	9667	Surgical Clipper Starter Kit	1	Case

## 3M™ Surgical Clippers with Pivoting Head and Clipper Blade Assembly

### Description

Rechargeable and corded surgical clippers for removing body hair from patients in preparation for surgery. Both models use disposable, single-use blade assemblies.



### Product Benefits

- April 1999 CDC Guidelines recommends that if hair removal is required, the use of a surgical clipper immediately before an operation is one means to help reduce the risk of surgical site infections.
- Effectively clips all types of body hair to eliminate need for multiple-style clippers for body hair removal.
- Rechargeable model provides approximately 20 minutes of run time from an 8-hour recharge, or 30-40 minutes of run time from a 12-14 hour recharge, allowing you to clip multiple patients from a single charge.
- Corded model provides continuous operation for longer preps and is a dependable backup for the rechargeable clipper.
- Single-use blades eliminate cross-contamination and are easy to snap on for use and pop off for disposal.
- Lightweight design fits hand comfortably, with blade at precise angle for easy clipping. Handle is easy to clean and disinfect to save time.

### Suggested Applications

- Removing body hair from patients in preparation for procedures where hair removal is desired


	Product Number	Product Information	Items/Box	Billing Unit
<b>Packaged Sterile</b>				
	9602	Surgical Clipper by Remington - rechargeable	1	Each
	9603	Surgical Clipper by Remington - cord operated	1	Each
	9600	Surgical Clipper by Remington Single-use Blades	40	Case

# Surgical Clippers and Prepping Solutions

## 3M™ DuraPrep™ Surgical Solutions

### Description

A patient skin prepping solution in a self-contained applicator that combines two proven broad-spectrum antimicrobials: alcohol for quick kill and fast drying, and iodophor for long-lasting antimicrobial protection.



### Product Benefits

- Achieves a 99.999% microbe kill (in-vitro) including antibiotic-resistant bacteria.
- Dries to a water-insoluble film that resists washing away by blood and saline.
- Self-contained, unit-dose applicator is designed for easy use and optimum coverage. Regulates flow of solution to minimize drips and mess.
- Simplified application allows for faster skin prep compared to traditional preps, to save you time.
- Maintains low bacterial counts under dressings for up to three days for long-lasting protection.
- Enhances drape adhesion to help assure that the sterile surface is not being compromised by incise drape lift.

### Suggested Applications

- The 8630 large-size applicator (26 ml) covers an area from neck to groin and is appropriate for procedures such as thoracic and CV, spine, total hip and knee, C-section, and major abdominal surgeries.
- The 8635 mid-size applicator (6 ml) covers at least an 8" x 10" area which is ideal for procedures requiring less coverage such as angiography, head and neck procedures, hand, foot and other minor procedures.

**DuraPrep solution contains alcohol which is flammable. To prevent fire, allow to dry completely (2-3 minutes on skin). Refer to package insert included with applicator for complete directions for use, warnings and precautions.**

	Product Number	Product Information	Applicators per Case	Billing Unit
--	----------------	---------------------	----------------------	--------------

#### Non-sterile - Bulk

	8630NS	Applicator (26 ml) 100 cotton-tipped swabs	50	Each
	8635NS	Applicator (6 ml)	100	Each

**Sterilization Note: May be EO processed 2 cycles. Do not gamma irradiate.**

	Product Number	Product Information	Applicators per Case	Billing Unit
--	----------------	---------------------	----------------------	--------------

#### Packaged Sterile

	8630	Applicator (26 ml)	20	Case
	8635	Applicator (6 ml)	50	Case

**Sterilization Note: May be EO processed. Do not gamma irradiate.**