

A History of 3M™ Tegaderm™ Brand



www.3M.com/Tegaderm



1978

3M begins work on inventing 3M™ Tegaderm™ Film Dressing

1981

Test launch in the U.S. First-ever with frame delivery system

1980

1982–1984
Tegaderm™ Film Dressing launched globally



1983

Optimized adhesive for universal applications: I.V. and wound

1983

First randomized, controlled clinical study

1985

1987

3M becomes U.S. market share leader in transparent film dressings

1988

Sales outside of U.S. surpass U.S. sales

1990

1986

First adhesive film dressing containing antimicrobials in adhesives

1991

3M introduces another, unique film and adhesive, Tegaderm™ HP Dressing for Holding Power in moist conditions



1993

3M invents unique dressing design with improved, simplified frame delivery system

1995



2003

Tegaderm™ Pediatric Dressing launched with printed teddy bears

2000

1996

3M sets a new standard in advanced film dressings, featuring borders, notches and enhanced securement




2005

2005

Tegaderm™ brand is extended to a family of trusted dressings for advanced wound care

2008
Launch of Tegaderm™ CHG Dressing



2010

2010

Tegaderm™ High Performance Foam Adhesive Dressing launched

2010

Global introduction of Tegaderm™ I.V. Advanced Securement Dressing, providing enhanced catheter securement and patient comfort

3M continues developing innovative products under the Tegaderm™ brand

2011

Tegaderm™ Diamond Pattern Dressing designed for use in moist environments

Celebrating 30 Years of Mutual Admiration

Every day, we are inspired by all you do to give your patients the best possible care. We value the trust you place in us and are committed to continually building that trust by providing people, products and ideas that raise the standard of care.

For over 30 years, 3M has collaborated with health care professionals around the world to develop products that simplify and improve care practices. That's how we developed Tegaderm™ transparent film dressings in 1982. Today, the Tegaderm™ brand offers an array of products for IV site and wound care, as well as the training, education and professional development resources you need, when and where you need them.

For more information about the Tegaderm™ brand family of dressings, visit www.3M.com/Tegaderm.





**Skin & Wound Care Division
3M Health Care**

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA

1-800-228-3957
www.3m.com/healthcare

3M and Tegaderm are trademarks of 3M.
Please recycle. Printed in U.S.A.
© 3M 2007, 2012. All rights reserved.
70-2009-9563-0