



BioTrust®


Finger or Face Computer Logon

3M Cogent's BioTrust® biometric logon software is a fast, easy, and secure way to replace your Windows® computer logon. Your face or fingerprint replaces the traditional username and password logon. Just sit down in front of your computer or swipe your finger across the fingerprint sensor and BioTrust will do the rest. Lighting conditions, such as side lighting, backlit areas, low ambient lighting, or lighting from the LCD screen itself are not an issue with 3M Cogent's advanced facial matching algorithm. All you need is a web camera or a fingerprint sensor.

The BioTrust console easily allows users to customize security level settings, enroll facial/fingerprint images, and manage Web SSO (Single Sign-On) subscriptions. Enrollment is a breeze with BioTrust's intuitive enrollment wizard: visual cues and videos guide you through the enrollment process. Your computer is biometrically secured in a matter of minutes!

BioTrust's automatic Web SSO functions as a password manager for website logons. Whenever you navigate to a secured Web site (Web mail, social networking sites, online banking) BioTrust automatically captures the logon information while keeping you in full control. The Dashboard with QuickLaunch allows quick and easy access to Web sites that have been enrolled in the Web SSO with one click, in one convenient location.

BioTrust is available in two configurations to meet your individual security needs: facial login and fingerprint or facial login.





BioTrust

Technical Specifications and Features

Features

Automatic Web SSO logon capture
Subscription Manager for SSO websites
Dashboard with QuickLaunch
SSO Scripting

Available Configurations

Face only
Fingerprint or Face

Supported Operating Systems

Windows® XP 32
Windows® Vista 32/64
Windows® 7 32/64

Minimum Hardware Requirements

Processor	Intel Atom leveled
RAM	1 GB
Hard Disk	100 MB free space
Camera	320 x 200
Fingerprint Sensor (not provided)	Validity Swipe

Recommended Hardware Requirements


Processor	Intel Centrino 1.6 GHz or higher
RAM	2 GB
Hard Disk	500 MB free space
Camera	Logitech with "RightLight"
Fingerprint Sensor (not provided)	Validity Swipe



► BioTrust's Automatic Web SSO allows the user to manage secure Web site logons.



► Easily navigate BioTrust's features and settings from the convenient console.



► Seven facial images are captured, providing greater accuracy.



► Intuitive, user-friendly interface guides the user through the enrollment process.