



3MSM CodeRyteSM CodeCompleteSM

- Leverages 3M's leading-edge natural language processing (NLP) technology to support physician documentation review and medical code assignment
- Helps eliminate coding overtime and backlogs
- Provides operational reports through an interactive dashboard that tracks productivity, volume, audit trails and turnaround times

The 3M advantage

3M CodeRyte CodeComplete is a full-service, outsourced coding solution comprised of a talented team of certified coders who can provide medical and surgical coding for both technical and professional services.

With 3M CodeRyte CodeComplete, organizations can have a 48-hour turnaround time for all of their medical records. In addition to improved coder productivity, clients can experience substantial savings through reductions in turnaround time, denials and days in A/R.

The challenge

Is your organization experiencing a shortage of strong, certified coders? Is your backlog out of control? Do you lack the time to train and manage new staff? Are you expecting short-term changes in staffing due to maternity leave or extended vacations?

The 3M solution: Certified coders at the ready

3M CodeRyte CodeComplete is a full-service, outsourced solution for the coding challenges your organization may face today or tomorrow. The strong roster of coders can help reduce coding backlog and improve coding turnaround time. Data entry is automated with outbound interfaces to billing systems, helping eliminate paper-based workflows and associated errors.

With **3MSM CodeRyteSM CodeAssistSM** as the foundational technology, the 3M CodeRyte CodeComplete team adheres to a coding-compliant workflow. The leading edge NLP technology in 3M CodeRyte CodeAssist supports physician documentation review and medical code assignment. NLP can be more consistent in its identification of medical facts than human coders and also help medical coders provide more consistent coding results.

3M CodeRyte CodeComplete clients receive both thorough retrospective and concurrent auditing capabilities. The flexible service also provides a seamless transition to build business without adding onsite staffing or handling the hassles of hiring, training and managing. 3M CodeRyte CodeComplete's per-report pricing model streamlines results, in lieu of more traditional hourly rates that hinge on productivity.



Results when you need them

With 3M CodeRyte CodeComplete, organizations can have a 48-hour turnaround time for all of their medical records. But improved coder productivity is just the start. 3M clients can experience substantial savings through reductions in turnaround time, denials and A/R days.

3M's talented team of certified coders provide medical and surgical coding for both technical and professional services. They deliver expert assignment of CPT[®], ICD and associated modifiers, DRGs, APCs, PQRS codes and more. With 3M CodeRyte CodeComplete, clients also receive regularly scheduled calls with internal coding experts to discuss coding findings and dictation issues for improved documentation.

Quality assurance

3M practices a comprehensive approach to quality assurance, which includes rigorous testing and screening of prospective coders, on-going auditing of coders, and training that is internal and external, generic and targeted.

3M's philosophy leverages a data-driven management style that facilitates consistent adherence to all relevant coding guidelines, client organization processes, and 3M's own determination of best practices based upon its work with well-respected organizations across the country. Additionally, our technology enables clients to do their own prospective and retrospective auditing when staffing permits.

Physician monitoring and operational reporting

3M CodeRyte CodeComplete clients have a new window into physician documentation. 3M's reporting tools easily identify areas for documentation improvements. At the same time, the service offers reporting on medical necessity edits by ordering/rendering physician, evaluation and management (E&M) element counting, bell curve service level reporting on a per-physician basis, and physician productivity by RVU.

3M also manages operational reports for 3M CodeRyte CodeComplete clients through an interactive dashboard that tracks productivity, volume, audit trails and turnaround time.

Consulting and ICD-10 preparation

3M CodeRyte CodeComplete clients have the option of including additional consulting services such as PQRS documentation and reporting, and physician trending compared to peer physicians. Organizations may also opt to include 3M's ICD-10 preparation services, which include:

- The ability to code ICD-9 and review parallel ICD-10 mapping during the auditing process
- Deliver customized physician training based on current documentation
- Perform automated mapping of ICD-9 or ICD-10 codes to bill by specific payer

Call today

For more information on how 3M solutions can assist your organization, contact your 3M sales representative, call us toll-free at **800-367-2447**, or visit us online at **www.3Mhis.com**.



Health Information Systems

575 West Murray Boulevard
Salt Lake City, UT 84123
U.S.A.
800 367 2447
www.3Mhis.com

3M, CodeRyte, CodeAssist and CodeComplete are service marks of 3M Company. CPT is a registered trademark of the American Medical Association.

Please recycle. Printed in U.S.A.
© 3M 2013. All rights reserved.
Published 1/13
70-2009-9213-2