

3M Food Safety
3M™ Electronic Pipettor II



Reliability & Quality Maximized

The 3M™ Electronic Pipettor II combines efficiency, versatility, and precision into one simple yet sophisticated tool. The 3M Electronic Pipettor II is pre-programmed to perform the most common 3M™ Petrifilm™ Plate dilutions and provides high accuracy and precision handling for micro-volumes of liquids. The microprocessor-based system reduces the possibility of human error and instrument contamination by controlling pipetting speed, calibration and volume selection. It's lightweight and ergonomic controls take the effort out of pipetting to help reduce the risk of repetitive strain injuries that are frequent with manual pipetting.

3M

Combine the cost-efficiency and labor-saving technology of 3M™ Petrifilm™ Plates with the features of the 3M Electronic Pipettor II

Enhanced productivity

- Pre-programmed to perform common dilutions within the tip
- Open programming mode for user-set pipetting and diluting
- Multiple plating with the multi-dispense feature on 5mL Electronic Pipettor II
- Auto Function Test at Start Up for selecting previous mode

Increased accuracy

- Microprocessor-controlled system
- Accurate, reproducible results
- Automatically detects calibration errors
- Nine speeds for aspirating and dispensing accuracy
- CE marked per EC Council Directive 2004/108/EC

Improved ergonomics

- Ergonomic design; reducing likelihood of repetitive strain injuries
- Operates either right- or left-handed
- Push-button auto tip ejector
- Soft-touch operating buttons
- Weighs less than 171 grams

3M recommends that users of all pipettors observe good ergonomic work practices including good workstation design and seating, as well as behavioral measures such as regular work breaks and varied work postures.

Ordering Information

Product	Cat. No.	Units/case
3M Electronic Pipettor II — 1mL	6502E	1
3M Electronic Pipettor II — 5mL	6503E	1
3M Pipettor Tips — 1mL, wide bore, racked, pre-sterilized	6486	960/case
3M Pipettor Tips — 1mL, wide bore, racked	6485	960/case
3M Pipettor Tips — 1mL, wide bore, bulk	6484	1000/case
3M Pipettor Tips — 5mL, wide bore, racked, pre-sterilized	6489	500/case
3M Pipettor Tips — 5mL, wide bore, bulk	6487	250/case



Important Information: The customer must determine whether this product is appropriate for use with its products and processes. Claims based on the absence of warranted characteristics or features can only be made, if in the particular case a concrete characteristic or feature has been warranted explicitly and in writing by the supplying 3M company. Besides this, all warranty and liability claims are subject to the Terms and Conditions of Sale of the supplying 3M company.



3M FOOD SAFETY
c/o 3M Deutschland GmbH
Carl-Schurz-Straße 1
41453 Neuss
Germany
Phone +49 (0) 2131/14-3000
Fax +49 (0) 2131/14-4443
Internet www.3m.com/micobiology

3M and Petrifilm are trademarks of 3M.
Please recycle. Printed in Germany
© 3M 2009. All rights reserved.
1427-101-EU 08/09