

3M™ Scotchcast™ One-Step Splints

Description

A splint that incorporates foam padding, splint and breathable backing for easy, one-step application. Prepackaged in a variety of sizes to meet application needs.

Product Benefits

- All-in-one splint eliminates the need for splint assembly for easy application and time savings. Faster application increases patient turnaround.
- Conformable and moldable for better patient fit.
- Provides strong, lightweight support to encourage early patient mobility.
- Available in a variety of pre-cut sizes to accommodate patient sizes and reduce the need to cut and customize splints.



Suggested Applications

- For use in the construction of common orthopedic splints

	Product Number	Size	Splints per Case	Billing Unit
Non-sterile				
	76210	2" x 10" (5 cm x 26 cm)	5	Case
	76312	3" x 12" (7,5 cm x 30 cm)	5	Case
	76415	4" x 15" (10 cm x 37,5 cm)	5	Case
	76430	4" x 30" (10 cm x 75 cm)	5	Case
	76530	5" x 30" (12,5 cm x 75 cm)	5	Case
	86335	3" x 35" (7,5 cm x 87,5 cm)	5	Case
	86545	5" x 45" (12,5 cm x 112,5 cm)	5	Case

3M™ Reinforcing Strips

Description

A reinforcing strip designed to aid in creating strong and durable fiberglass casts, particularly those likely to be exposed to significant challenges of breakdown or stress.

Product Benefits

- Reinforcing strips are easy to use for fast reinforcement of high stress casts.
- Provide strength and stability to improve cast durability and wear.
- Wide weave fiberglass mesh allows the cast to breathe better than with closed weave splints for cooler, lighter casts.



Suggested Applications


- Reinforcing fiberglass casts, especially those likely to be exposed to significant challenges of breakdown or stress

	Product Number	Size	Strips per Case	Billing Unit
Non-sterile				
	70941	4" x 18" (10 cm x 46 cm)	10	Case
	70943	4" x 36" (10 cm x 92 cm)	10	Case

3M™ Synthetic Cast Paddings

Description

A non-woven polyester padding that is porous, highly conformable, and strong yet easy to tear. Suitable for use with all casting materials.



Product Benefits

- Spun Dacron polyester padding provides protection and patient comfort, while helping reduce water retention for a drier cast.
- Strong and durable padding is easy to tear.
- Conforms easily to body contours.
- Porous to allow the skin to breathe.
- Versatile padding can be used with all casting materials, plaster and synthetic.
- Available in a variety of widths to meet your casting needs.

Suggested Applications

- Padding under casts

	Product Number	Size	Rolls per Case	Billing Unit
Non-sterile				
	CMW02	2" x 4 yds (5 cm x 3,6 m)	80	Case
	CMW03	3" x 4 yds (7,5 cm x 3,6 m)	80	Case
	CMW04	4" x 4 yds (10 cm x 3,6 m)	80	Case
	CMW06	6" x 4 yds (15 cm x 3,6 m)	80	Case

3M™ Synthetic Cast Stockinets

Description

An ideal component of a high-quality cast, this synthetic stockinet is intended for use as an underlayment for all standard casting applications. Suitable for use with all casting materials, plaster and synthetic.



Product Benefits

- Stockinet helps reduce water retention and improves patient comfort.
- 100% knitted polyester stockinet provides conformability for easy use.
- Can be used for all standard casting applications and with all casting materials for greater versatility.
- Available in a variety of widths to meet your casting needs.

Suggested Applications

- For use as an underlayment for all standard casting applications

	Product Number	Size	Rolls per Case	Billing Unit
Non-sterile				
	MS01	1" x 25 yds (2,5 cm x 22,8 m)	1	Case
	MS02	2" x 25 yds (5 cm x 22,8 m)	1	Case
	MS03	3" x 25 yds (7,6 cm x 22,8 m)	1	Case
	MS04	4" x 25 yds (10,1 cm x 22,8 m)	1	Case
	MS06	6" x 25 yds (15,2 cm x 22,8 m)	1	Case
	MS08	8" x 25 yds (20,3 cm x 22,8 m)	1	Case