

3M™ Red Dot™ Diaphoretic Foam Monitoring Electrodes

Description

A rectangular-shaped electrode that offers a high-tack adhesive and foam backing designed for diaphoretic clinical situations.



Product Benefits

- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- Use of built-in abrader disc improves trace quality for better results.
- 3M™ Moistac™ Adhesive provides greater adhesion to the skin, especially in moist conditions, which results in better trace quality.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.
- Foam backing is occlusive to fluids and is easy to handle and remove.

Suggested Applications

- General ECG monitoring situations where skin moisture is an issue
- For use in ER, OR, Stress Test, Cardiac Rehab and Echo-Cardiography departments

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2230	Foam tape backing, adult size, diaphoretic	50	20	1,000	Bag
	2270	Foam tape backing, adult size, diaphoretic	25	40	1,000	Bag

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2270-3	Foam tape backing, adult size, diaphoretic, 3 electrodes/bag	20	10	600	Case
	2270-5	Foam tape backing, adult size, diaphoretic, 5 electrodes/bag	20	10	1,000	Case

3M™ Red Dot™ Diaphoretic Soft Cloth Monitoring Electrodes

Description

A rectangular-shaped electrode that offers a high-tack adhesive and soft cloth backing designed for diaphoretic clinical situations.



Product Benefits

- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- Use of built-in abrader disc improves trace quality for better results.
- 3M™ Moistac™ Adhesive provides greater adhesion to the skin, especially in moist conditions, which results in better trace quality.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.
- Permeable soft cloth backing releases perspiration, making it ideal for long-term monitoring in diaphoretic situations.

Suggested Applications

- General ECG monitoring situations where skin moisture is an issue
- For use in ICU/CCU, Stress Telemetry, and Holter Monitoring departments

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2231	Soft cloth tape backing, adult size, diaphoretic	50	20	1,000	Bag
	2271	Soft cloth tape backing, adult size, diaphoretic	25	40	1,000	Bag

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2271-3	Soft cloth tape backing, adult size, diaphoretic, 3 electrodes/bag	20	10	600	Case
	2271-5	Soft cloth tape backing, adult size, diaphoretic, 5 electrodes/bag	20	10	1,000	Case

3M™ Red Dot™ Foam Monitoring Electrodes

Description

Tear-drop style monitoring electrodes with foam backing.

Product Benefits

- Foam backing is comfortable for patients, occlusive to fluids, and easy to handle for quick application.
- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.
- Use of built-in abrader disc (2259) improves trace quality for better results.



Suggested Applications

- ECG monitoring
- For use in ER, OR, Stress Test, Cardiac Rehab, and Echo-Cardio departments

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	---------------	---------------------	--------------

Non-sterile

	2237	Foam tape backing, adult size	50	20	1,000	Bag
	2259	Foam tape backing, adult size, abrader	25	40	1,000	Bag

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	----------------	---------------------	--------------

Non-sterile

	2237-3	Foam tape backing, 3 electrodes/bag	20	10	600	Case
	2237-5	Foam tape backing, 5 electrodes/bag	20	10	1,000	Case
	2259-3	Foam tape backing, abrader, 3 electrodes/bag	20	10	600	Case

3M™ Red Dot™ Monitoring Electrodes with Foam Tape and Sticky Gel



Description

Multi-purpose monitoring electrodes with foam backing and sticky gel are ideally suited for the OR, ER and diagnostic stress testing applications.

Product Benefits

- Conveniently packaged in individual strips of ten electrodes for easy use and less waste.
- Remains fresh for 30 days out of the bag.
- High performance adhesive offers instant adhesion and up to five-day wear with excellent patient skin care and comfort.
- Sticky gel delivers excellent trace quality resulting in better patient information.
- Smaller size electrode is easier to place, and minimizes skin irritation for better patient comfort.
- Microreplicated Trace Prep surface (2570) gently prepares skin for optimal trace quality.
- 2570 electrode is radiolucent to allow for easy, cost-effective use throughout the hospital.

Suggested Applications

- ECG monitoring
- For applications requiring high adhesion
- For hospital-wide use; especially suited for the ER, OR and Catheter Lab or ICU

	Product Number	Product Information	Strips per Bag	Bags per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2560	Foam tape backing, sticky gel, 10 electrodes/strip	5	20	1,000	Case
	2570	Foam tape backing with radiolucent stud and trace prep, 10 electrodes/strip	5	20	1,000	Case
	2570-3	Foam tape backing with radiolucent stud and trace prep, 3 electrodes/strip	1	200	600	Case
	2570-5	Foam tape backing with radiolucent stud and trace prep, 5 electrodes/strip	1	200	1,000	Case

3M™ Red Dot™ Monitoring Electrodes with Micropore™ Tape Backing



Description

A round-style monitoring electrode with 3M™ Micropore™ Surgical Tape backing.

Product Benefits

- Micropore backing is breathable, flexible and conformable, yet is gentle to the skin for increased patient comfort.
- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.
- Use of built-in abrader disc improves trace quality for better results (only on 2249).
- Radiolucent to minimize the need to replace or remove electrodes in most radiological imaging applications.

Suggested Applications

- ECG monitoring situations that require a comfortable, breathable and permeable backing
- For use in ICU/CCU, Telemetry, and Holter Monitoring departments

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	---------------	---------------------	--------------

Non-sterile

	2239	Micropore tape backing, adult size	50	20	1,000	Bag
	2249	Micropore tape backing, adult size, abrader	25	40	1,000	Bag

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	----------------	---------------------	--------------

Non-sterile

	2249-3	Micropore tape backing, abrader, 3 electrodes/bag	20	10	600	Case
--	--------	---	----	----	-----	------

3M™ Red Dot™ Pediatric Monitoring Electrodes with Micropore™ Tape Backing



Description

A traditional-style monitoring electrode with 3M™ Micropore™ Surgical Tape backing that is ideal for pediatric applications.

Product Benefits

- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- Micropore tape backing is breathable, stretchable and conformable, yet is gentle to the skin for increased patient comfort.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.

Suggested Applications

- For use in pediatric units, including Stepdown/ Observation and Pediatrics PICU

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2248	Micropore tape backing, pediatric size	25	40	1,000	Bag

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2248-3	Micropore tape backing, pediatric size, 3 electrodes/bag	20	10	600	Case

3M™ Red Dot™ Soft Cloth Monitoring Electrodes

Description

A round-style monitoring electrode with soft cloth backing for general ECG monitoring.

Product Benefits

- Soft cloth backing is breathable, flexible and conformable, yet is gentle to the skin for increased patient comfort.
- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.
- Use of built-in abrader disc improves trace quality for better results (only on 2255).



Suggested Applications

- ECG monitoring
- For use in ICU/CCU, Telemetry, and Holter Monitoring departments

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	---------------	---------------------	--------------

Non-sterile

	2238	Soft cloth tape backing, adult size	50	20	1,000	Bag
	2255	Soft cloth tape backing, adult size, abrader	25	40	1,000	Bag

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	----------------	---------------------	--------------

Non-sterile

	2238-3	Soft cloth tape backing, 3 electrodes/bag	20	10	600	Case
	2238-5	Soft cloth tape backing, 5 electrodes/bag	20	10	1,000	Case
	2255-3	Soft cloth tape backing, abrader, 3 electrodes/bag	20	10	600	Case

3M™ Red Dot™ Infant Soft Cloth Monitoring Electrodes



Description

A traditional-style monitoring electrode with soft cloth backing for infant applications.

Product Benefits

- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- Soft cloth backing is breathable, stretchable and conformable, yet is gentle to the skin for increased patient comfort.

Suggested Applications

- For use in neonatal units including Stepdown/Observation

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2258-3	Soft cloth tape backing, infant size, single pack, 3 electrodes/bag	20	10	600	Case

3M™ Red Dot™ Clear Plastic Monitoring Electrodes

Description

A tear-drop style monitoring electrode with clear plastic backing that is ideal for situations when immediate skin assessment is required without the need to remove the electrode.



Product Benefits

- Clear plastic backing allows for easy viewing of patient's skin throughout wear cycle which minimizes the need to remove or replace electrodes when patients move between departments.
- Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces.
- The gel's low chloride content minimizes skin irritation for increased patient comfort.
- 3M™ SureSeal™ Cap Seal prevents premature product dryout to help ensure product freshness for 30 days out of the bag.

Suggested Applications

- ECG monitoring
- For use in ER, OR, Stress Test, Cardiac Rehab, and Echo-Cardio departments

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	---------------	---------------------	--------------

Non-sterile

	2235	Clear tape backing, adult size	50	20	1,000	Bag
--	------	--------------------------------	----	----	-------	-----

	Product Number	Product Information	Electrodes per Box	Boxes per Case	Electrodes per Case	Billing Unit
--	----------------	---------------------	--------------------	----------------	---------------------	--------------

Non-sterile

	2235-3	Clear tape backing, 3 electrodes/bag	20	10	600	Case
	2235-5	Clear tape backing, 5 electrodes/bag	20	10	1,000	Case

3M™ Red Dot™ Resting EKG Electrodes



Description

Tab style resting electrodes for 12-lead EKG applications.

Product Benefits

- Tab-style electrode has pressure-sensitive conductive adhesive that is easy to apply and remove, and provides rapid and intimate skin contact to ensure high quality trace pickup.
- Radiolucent to minimize the need to replace or remove electrodes in most radiological imaging applications.
- Extra large tabs for easy handling and attachment.
- Semi-stiff backing reduces roll-off and minimizes electrode lift.

Suggested Applications

- 12-lead EKG resting diagnostic monitoring

	Product Number	Product Information	Electrodes per Bag	Bags per Case	Electrodes per Case	Billing Unit
Non-sterile						
	2360	Resting tab electrode 10 electrodes/card, 10 cards/bag	100	40	4,000	Case